

MAP INDUSTRIES has established itself nationally as a pioneer in innovative Equipments and Solution with Futuristic vision in hospital CSSD.

HPHV FULLY AUTOMATIC AUTOCLAVE STERILIZER.

Sterilization at standard 121°C and 134°C cycles. | On-screen display of phase current & set values, temperature, pressure & humidity. | Vacuum pulsing and drying to achieve low Sterility Assurance Level (SAL). | Automatic control of temperature, moisture, values and all critical parameters. | Reports auto generated in pen drive, hard disk or printer.

HPHV Rectangular Autoclave



Rectangular on Stand

- Gage
- Compressed Air Not Required
- HMI
- Printer & Emergency Stop
- Temperature & Humidity Sensor
- USB Port
- SS 316L



Sliding Door Compact Modal



Hinge Door Compact Modal

S.No.	Volume (Cuft)	Liter
1.	16.00	450
2.	20.00	550
3.	31.00	850

S.No.	Dimension (in)	Liter	Drum
1.	16 Ø x 24 L	81	2
2.	20 Ø x 36 L	185	3
3.	20 Ø x 48 L	248	4

HPHV Cylindrical Autoclave



Cylindrical on Stand



Cylindrical Compact Modal

- Gage
- Compressed Air Not Required
- HMI
- Printer & Emergency Stop
- Temperature & Humidity Sensor
- USB Port
- SS 316L



Cylindrical Vertical Autoclave

B type, Flash Sterilizer

- » B type fully automatic front loading autoclave with a capacity of 22 Liter.
- » With vacuum, dry cycle & micro processor controlled.
- » Double safety lock on the door.
- » Can be used at 121 C / 134 C during sterilization.
- » 2 Storage tanks for fresh & used water.
- » Wrapped/Unwrapped, Porous, Solid Instruments including Cotton Swabs can be sterilized.
- » Bowie & dick test qualified, Inbuilt sensor for distilled water.
- » Digital display of pressure and temperature.
- » Chamber Size: 250 mmx 450 mm



Flash Autoclave

Ultrasonic Washer



Ultrasonoc Washer

- » ULTRASONIC CLEANER For Surgical Instruments
- » Micro Controller base Digital Timer .
- » With latest MOSFET base technology.
- » Automatic frequency tuning for maximum output.
- » Low transducer losses permit maximum cleaning action.
- » Modular generator for easier modification/servicing.
- » Rugged and easy to clean stainless steel housing.
- » Extensively protected electronics circuits for longer and safe operation.
- » Frequency 33 OR 40KHz + 3 KHz
- » Input voltage range of 170AC-270AC.50 Hz
- » Inbuilt Digital Timer
- » Auto Degas

Other Products



Sealing Machine



Fogger Machine



All type of indicators

MAP INDUSTRIES

CHANNEL PARTNER

CSSD TURNKEY PROJECT

- 📍 Our Address: A-103,104, Swarnim Industrial Park, Dhamatvan Road, Bujrang Village, Bakrol, Ahmedabad, Gujarat - 382433, India.
- 🌐 Website: www.mapindus.com
- ✉ Email: info@mapindus.com, mapindustries94@gmail.com
- ☎ Abdul Hanif: 9687341157 | 9687771157