

**MAP INDUSTRIES** has established itself nationally as a pioneer in innovative Equipments and Solution with Futuristic vision in hospital CSSD.

## Fully Automatic ETO (Ethylene Oxide) Sterilizer

on-screen display (FUJI) of phase, current & set values, temperature, pressure & humidity.  
automatic control of temperature, moisture, vacuum & valves.  
current/previous reports printed through printer.  
PLC (Honey well) based automatic with 3 (Cold, Hot, Medium) program cycles  
Material using in chamber SS material 316L.  
Auto/manual operation are available , Negative pressure cycle.

## Fully Automatic PLC Based ETO Sterilizer

S.No.	Volume (Cuft)	Liter	ETO Cartridge (gm)
1.	2 cuft	55	40
2.	4 cuft	110	100
3.	8 cuft	220	170




Polymex Tpe



Medical Instruments



Long Lumens



Surgic Instruments



Medical Instruments



Implants Items



Fully Automatic Full Body Model



Fully Automatic Table Top Model



Sealing Machine



ETO Cartridge



All Type Indicators



Sterilization Roll

# Central Supply & Sterilization Department (CSSD) Furniture



**SS ENCLOSED DISTRIBUTION TROLLEY FOR STERILE GOODS**



**Ss Washing Sink**



**SS PACKING, SEALING & LABELING**



**LAPAROSCOPY TROLLEY**



**SS Table Trolley**



**SS Linen Inspection Table**



**SS STOOL**



**Storing Rack**



**SS FOOT OPERATED SINK - 2 BAY**



**SS FOOT OPERATED SINK - 1 BAY**



**SS Modular Basket**



**SS Scope Tray**



**SS Pass Box**



**SS Dressing Drum**



**SS Trolley**



**SS Crash Cart Trolley**



**SS Water Gun**



**SS Perforated Instrument Box**



**SS wall mounted Storing Rack**



**Three Tyre Trolley**

**MAP INDUSTRIES**

**CSSD TURNKEY PROJECT**

📍 Our Address: A-103,104, Swarnim Industrial Park, Dhamatvan Road, Bujrang Village, Bakrol, Ahmedabad, Gujarat - 382433, India.

🌐 Website: [www.mapindus.com](http://www.mapindus.com)

✉ Email: [info@mapindus.com](mailto:info@mapindus.com), [mapindustries94@gmail.com](mailto:mapindustries94@gmail.com)

☎ Abdul Hanif: 9687341157 | 9687771157

**CHANNEL PARTNER**